# Home Care Technical Specification Tool Free Transfer Bench DN7065


## **Stock Numbers:**

DN7065 – Tool-Free Transfer
 Bench
 White

#### **Materials:**

**Tool-Free Transfer Bench** is made of ABS Plastic, Aluminum, EPDM, and Stainless.

## Dimensions:



### **Installation Instructions:**

- Determine orientation of chair relative to water supply in order for transfer bench to be positioned outside of bathtub.
- Locate the two under seat support beams: Hold the pieces by the straight end with the curved end away from you.
- 3. Locate the flat, indented area of each beam: Place the two flat areas together, so they are touching, and the two beams form an "X" shape.
- 4. Locate the leg unit with top posts that point inward: Align the posts with the straigh ends of the crossed support beams on the small end of the "X".
- 5. Slide the top posts of the leg unit over the ends of the crossed support beams.
- 6. Take the remaining leg unit, with the top posts that point straight out, and align it with the curved end of the support beams, at the large side of the "X".
- 7. Slide the top posts of the leg unit over the ends of the crossed support beams.
- 8. Locate the transfer bench. Place it seat side down.
- Place the curved ends of the crossed support beams into the grooves on the underside of the transfer bench, aligning the holes in the support beams with the holes inside the grooves.
- Use the assembly knobs to attach the transfer bench to the support beams.
   Thread the washers onto the bolt, small lock washer first, followed by the cup washer.
- 11. Locate the large seat bottom. Place it under the support beams so that the center
- 12. of the "X" slides into the center of the "X" shaped grooves on the large seat bottom.

Stock Number:	TF Transfer Bench
By Creative Specialties	

25300 AL MOEN DRIVE • NORTH OLMSTED, OH 44070 CORPORATE HEADQUARTERS: (800) 321-8809 • FAX: (800) 848-6636 INTERNATIONAL: (440) 962-2000 • FAX: (440) 962-2726 CUSTOMER SERVICE (USA): 800-882-0116 • FAX: (888) 379-2720 WWW.HOMECARE.MOEN.COM

